

Play **BINGO** with
Dentsply Sirona



Road Map to SiroLaser Blue Demo



Probing Questions

Do you own an electrosurge?

Yes No

Do you own a laser unit?

Yes No

If yes, what brand?

Are you still using a scalpel for
gingival contouring?

Yes No

Do you do a lot of frenectomies and
gingivectomies in your clinic?

Yes No

Are you sending out any cases
that would need gingival
restructuring or soft tissue cutting?

Yes No

Why a Demo?

- ▶ SiroLaser Blue is equipped with a high-tech laser module, which provides 3 wavelengths in 1 device. This enables it to cover more than 20 indications.
- ▶ **Blue laser light (445 nm)** has a much higher absorption rate in soft-tissue (i.e. hemoglobin and melanin) than conventional infrared diode lasers. This leads to a much improved soft-tissue cutting efficiency which allows non-contact cutting. In a **CAD/CAM** workflow for example, it minimizes damage and bleeding to the tissue, resulting in clearly defined margins around the preparation site. Therefore, it eliminates the use of retraction cords.
- ▶ **Infrared laser light (970 nm)** is effective in germ and bacterial reduction for the **hygiene** and **perio** workflows. In an **endo** treatment, the use of infrared after the conventional treatment of the root canal leads to a better long-term prognosis and a reduced use of antibiotics.
- ▶ **Red laser light (660nm)** is used for **photobiomodulation** to improve wound healing and tissue regeneration, supporting a wide range of indications in both surgical and therapeutic procedures. It can be used in **TMJD treatment** for example, to reduce pain and improve mandibular movement.

TO BOOK A SIROLASER BLUE DEMO,

please complete the form below and submit it to the BINGO online platform.

Henry Schein Rep name:

Practice name:

Dentist:

Address:

Email:

Phone: