



# PracticeProtect™

The importance of protecting patients from exposure to infectious diseases is key to a safe dental visit, not to mention a lifetime of oral health. Throughout the course of your day, our PracticeProtect™ range of infection prevention products cover all the touch points that could potentially spread disease. It's a simple and highly effective system to protect everyone's safety—you, your patients and your practice.



## FlashTips®

### Disposable Air/Water Syringe Tips

- Rounded tip design helps prevent O-ring leakage or air/water contamination issues
- Locking indicator to help ensure proper tip placement
- Disposable, single-use tips minimize risk of cross-contamination



SKU	Description	Product Qty	Quantity
23071	FlashTips Disposable Air/Water Syringe Tips	250	
23075	FlashTips Disposable Air/Water Syringe Tips	1200	



## Purevac® SC


### Evacuation System Cleaner

Formulated specifically for systems with amalgam separators.

- Neutral pH formula meets guidelines for evacuation line cleaners
- Walkabout dispenser creates a turbulent mixture of air and solution to clean evacuation system with ease.



SKU	Description	Product Qty	Quantity
21132	Purevac SC Evacuation System Cleaner	2L	
21135	Purevac SC Evacuation System Cleaner	5L	
21116	Purevac SC Evacuation System Cleaner Dispenser	1	



## Resurge®

### Instrument Cleaning Solution

- Neutral pH, Biodegradable and environmentally friendly
- Removes difficult protein soil faster than leading instrument cleaning solution.
- Unique formula helps protect your instruments passive layer (outer coating) reducing the potential for stains, spots and even rust.



SKU	Description	Product Qty	Quantity
21520	Resurge® 24 packets	24 x 0.5 oz	
21521	Resurge 1L Bottle	338 oz	

## Purevac® HVE System

A single instrument for evacuating aerosols, splatter, fluid and debris.

- HVE removes 90% more aerosols generated during ultrasonic scaling compared to a low-volume saliva ejector.
- 135% greater removal of fluids than a low-volume saliva ejector\*\*
- Provides visibility, illumination and suction



SKU	Description	Product Qty	Quantity
11000	Purevac® HVE System Kit	1	
11010	Purevac® HVE Mirror Tips	3 pk	

For more information on PracticeProtect™ solution, visit <http://www.dentsplysirona.com/practiceprotect>

## Assure Plus® Sterilization Pouches



Self Seal Sterilization Pouches with unique peel-off internal monitor. Internal monitor is inserted inside pouch to confirm heat penetration according to CDC and IPAC guidelines. No need to order secondary strips.

SKU	Description	Dimensions	Product Qty	Quantity
83000	3.5" x 6.5"	(89mm x 165mm)	200	
83001	2.25" x 5"	(57mm x 127mm)	200	
83003	2.75" x 10"	(70mm x 254mm)	200	
83005	3.5" x 10"	(89mm x 254mm)	200	
83010	5.25" x 11"	(133mm x 279mm)	200	
83015	7.5" x 14"	(190mm x 355mm)	200	
83020	10" x 16"	(254mm x 406mm)	200	
83025	12" x 18"	(305mm x 457mm)	200	
83045	5.5" x 14"	(140mm x 355mm)	200	
83075	5.25" x 7.5"	(133mm x 191mm)	200	
83080	11" x 16"	(279mm x 406mm)	200	

## Disposa-Shield® Surface Disposable Barrier



- Disposable covers for high-touch dental equipment and surfaces
- Designed to help prevent cross-contamination

SKU	Description	Dimensions	Product Qty	Quantity
A88005	For Contra Angle Handpiece	8.9" Long	500	
A88006	For Low-Speed, Long Handpiece	7.6" x 1.5"	500	
A88007	For Air/Water Syringe	5.6" x 2.3"	500	
A88008	For Light Handle	6.8" x 2.95" x 1.4"	500	
A880095	For Hoses and Tubing	24" x 2"	175	
A880105	For Chair and X-Ray Head	15" x 27.5" (Expands to 27.5" x 18")	150	
A880115	Multi-Purpose Use Sheets	4" x 6"	Roll of 1200 Perforated Sheets	
A88013	For Etchant or Sealant Syringes	7.6" x 1.375"	500	
A88014	For Tray, Size 10.5" x 14.5"	10.5 x 14.5"	500	
A88015	For Tray, Size 12" x 17"	12" x 17"	500	
A88016	For Chair Headrest	13" x 13"	250	
A88019	For Cordless Prophy System Inner Module	Inner Module	500	



Don't forget to follow us on our Dentsply Sirona Preventive Youtube channel.

For more information about Dentsply Sirona, please visit [www.dentsplysirona.com/en-ca](http://www.dentsplysirona.com/en-ca) or contact your Dentsply Sirona Representative.

Name \_\_\_\_\_ Phone \_\_\_\_\_ Email \_\_\_\_\_